



## Disposable PE Laminated Shoe Cover

FDA CE MADE  
IN TURKEY

### Product Details

**Documents:** CE, ISO 13485

**Brand Name:** Dexxon Medical

**Model Nr:** DXNMD-DNSSC-02

**Certificate Nr:** NVA-EC-21022302

**Product Name:** Disposable PE Laminated shoe cover

**Use:** Disposable surgical shoe cover

**Material:** % 100 Spunbond PP non-woven shoe cover

**Material quality:** PE (polyethylene) film-coated non-woven fabric

**Fabric Type:** SMS fabric microporous polypropylene + polyethylene

**Usage:** Standard protective shoe cover

**Sterilization:** Non-sterile

**Weight:** 55 g / m<sup>2</sup>

**Size:** Standard size

**Colour:** Medical blue, white



## Disposable PE Laminated Shoe Cover



**Performance indicators:** Technical requirements for disposable laminated shoe cover  
EN 14126 :2003 ( Type 5-B, Type 6-B )

Tested in accordance with PPE standard Cat III ; Type 5/6

Type 5/6 Microscopus protective coverall with hood

SMS polypropylene base material with PE microporous film coating

**Related Directive :** Personal Protective Equipment 2016/25 (EU)

**Product features:** Non-woven non-sterile disposable medical protective coveralls

The product's resistance to synthetic blood penetration, breaking strength, and moisture permeability are above the national standard and have excellent performance for both safety and comfort.

Disposable and Breathable

Spunbond (SS or SMS) consists of 100% PP laminated non-woven fabric

Antistatic

Chemical resistant

Microporous polyethylene laminate

Comfort and durability guaranteed

**Performance indicators:** Technical requirements for disposable laminated shoe cover

**Application:** Hospital, medical, surgical, clean room,

**Precautions:**

The product should be stored in a clean, dry place with proper ventilation and out of direct sunlight at 0°C-30°C (its relative temperature should not exceed 75%)