



Disposable PP Shoe Cover

FDA CE MADE
IN TURKEY

Product Details

Documents: CE, ISO 13485

Brand Name: Dexxon Medical

Model Nr: DXNMD-DPPSC-03

Certificate Nr:

Product Name: Disposable PP shoe cover

Use: Disposable shoe cover

Material: % 100 Spunbond PP non-woven shoe cover

Material quality: PP non-woven fabric

Fabric Type: SS or SMS fabric

Usage: Standard protective coverall

Sterilization: Non-sterile

Weight:

Classification:

Certifications :

Size: Free Size

Colour: Blue



Disposable PP Shoe Cover

FDA CE MADE
IN TURKEY

Product features:

Made of polyethylene or polypropylene non-woven 30 g/m² or 40 g/m²
(for shoe covers with anti-skid layer)

Elastic rubber band at the opening

Protection against contaminating rooms where higher cleanliness

Standards are required

Increased durability

Disposable

Spunbond (SS or SMS) consists of 100% PP non-woven fabric

Antistatic

Dust-proof

Performance indicators: Technical requirements for disposable laminated
shoe cover

Application: Hospital, medical

Precautions:

The product should be stored in a clean, dry place with proper ventilation
and out of direct sunlight at 0°C-30°C

(its relative temperature should not exceed 75%).