



Disposable Surgical Gown SMS PE Laminated Reinforced, Sterile



Product Details

Product Description : SMS PE Laminated Reinforced Surgical Gown

Brand Name: Dexxon Medical

Model Nr: DXNMD- DNSLRSG-09

Certificated : CE and ISO13485

Certificate Nr:

Product Name: Disposable SMS Laminated Reinforced Sterile Surgical Gown

Material quality: PE (polyethylene) film-coated non-woven fabric

Type : Laminated Reinforced Surgical Gown Level 3 AAMI

Usage: Reinforced Surgical Gown

Sterilization: Sterile by ethylene oxide

Weight: 45 -55 g / m²

Classification: Class I

Regular Color: Medical Blue

Basic Info

Name: SMS PE Laminated Reinforced Surgical Gown

Color : Blue

Size : S.M.L.XL.XXL.XXXL

Regular Material: PP/SMMS



Disposable Surgical Gown SMS PE Laminated Reinforced, Sterile

  MADE
IN TURKEY

Performance indicators: Technical requirements for disposable polyethylene laminated reinforced isolation clothing products / EN 13795-1:2019; AAMI PB70: 2012

Level : Level 3

Feature:

SMS Non-Woven Fabric
Excellent water repellence to protect surgeons from contamination.
Prevention of cross contamination by reducing lint
Resistant to liquids with its barrier feature
Resistant to moisture and bacterial penetration
Puncture, tear and abrasion resistant.
Free from toxic substances
EO-Sterilized

Application:

This gown can be used during the surgery operations.
This product is disposable and it prevents the bacterial transfer between the medical team and the patient, reduces the occurring of possible contamination and infection between the patient and the medical team.

Precautions:

The product should be stored in a clean, dry place with proper ventilation and out of direct sunlight at 0°C-30°C (its relative temperature should not exceed 75%).



SMMS Surgical Gown, Disposable Patient Use SMMS Scrub Suit, Non-Sterile

FDA CE MADE
IN TURKEY

Product Details

Product Description : SMMS Surgical Gown / Scrub Suits

Brand Name: Dexxon Medical

Model Nr: DXNMD-DNSSS14

Certificated : CE and ISO13485

Certificate Nr: NVA-EC-21061104

Product Name: Disposable non-sterile short-sleeved patient gown

Type : SMS Surgical Gown / Scrub Suits

Usage: Standard patient gown

Sterilization: Non-sterile

Weight: 35-40 g / m²

Classification: Class I

Regular Color: Medical Blue

Basic Info

Name: SMS Surgical Gown / Scrub Suits

Color : Blue

Size : S.M.L.XL.XXL.XXXL

Regular Material: PP/SMMS