



## Disposable Surgical Gown SS, Sterile

FDA CE MADE  
IN TURKEY

### Product Details

**Product Description :** SS Surgical Gown

**Brand Name:** Dexxon Medical

**Model Nr:** DXNMD-DNSMG-SS04

**Certificated :** CE and ISO13485

**Certificate Nr:** NVA-EC-21052601

**Product Name:** Disposable sterile surgical gown

**Type :** SS Surgical Gown

**Usage:** Standard surgical gown

**Sterilization:** Sterile by ethylene oxide

**Weight:** 35-45 g / m<sup>2</sup>

**Classification:** Class IS

**Regular Color:** Medical Blue

### Basic Info

**Name:** SS Surgical Gown

**Color :** Blue

**Size :** S.M.L.XL.XXL.XXXL

**Regular Material:** PP/SS



## Disposable Surgical Gown SS, Sterile

FDA CE MADE  
IN TURKEY

### Performance indicators:

Technical requirements for disposable surgical gown  
EN 13795 : 2019

### Feature:

Disposable  
SS Fabric  
EO-Sterilized

### Application:

Our surgical gown is designed to be worn by healthcare professionals during surgical procedures.

Fabric feature is blue SS fabric.

The apron features a loop neck closure for an easily adjustable, gapless custom fit.

Cuffs at the wrists for added protection

### Precautions:

The product should be stored in a clean, dry place with proper ventilation and out of direct sunlight at 0°C-30°C (its relative temperature should not exceed 75%).