



**Disposable Surgical
Isolation Gown
SMS PE Laminated
Type 6-B Non-Sterile**

FDA **CE** **MADE
IN TURKEY**

Product Details

Product Description : SMS PE Laminated Surgical Isolation Gown

Brand Name: Dexxon Medical

Model Nr: DXNMD-DNSLIG13

Certificated : CE and ISO13485

Certificate Nr: NVA-EC- 21061103

Product Name: Disposable Isolation Clothing (Coveralls)

Type : Type 6-B Offering limited protective performance against liquid chemicals

Usage: Surgical Isolation Gown

Sterilization: Non-Sterile

Weight: 45 -55 g / m²

Classification: CAT III

Regular Color: Medical Blue

Basic Info

Name: SMS PE Laminated Surgical Isolation Gown

Color : Blue

Size : S.M.L.XL.XXL.XXXL

Regular Material: PP/SS or SMS



Disposable Surgical Isolation Gown SMS PE Laminated Type 6-B Non-Sterile



Performance indicators: Technical requirements for disposable isolation clothing polyethylene laminated extra reinforced products

EN 14126:2003 AC:2004 , EN 13034:2005 +A1:2009, EN 13795-1:2019; AAMI PB70: 2012

Type : Type 6-B

Level: Level 3

Feature:

Disposable laminated surgical isolation gown

Scope of Application: It is used for general isolation in clinics, wards, and laboratories of medical institutions. It is used for routine protection, laboratory protection, work and production resumption protection, and protection during public place disinfection

Structural Composition: The product is mainly made of non-woven cloth through cutting and sewing. It is provided in a non-sterile condition and only for single use.

Clothing Features and Overview:

The isolation gown is one-piece, coming in both medium and long styles, as shown in the figure. Rationally designed for free movement, it is easy to put on and take off. The sleeve cuffs have an elastic closure. It is dry, clean, mildew-free, and free of adhesions, cracks and cavities.

Humanized design: The waist ties with a knot, which is convenient for operating personnel of similar height but different waist size.

Comfortable fabric: New materials, soft, light, breathable.

Waist tie: Use of a knot-tie design allows for flexible adjustment, enhancing the clothing's personalized fit.

Cuff lock-stitching: Elastic cuff design ensures both heightened protection and free movement.

Characteristics

Excellent water repellence to protect surgeons from contamination.

Prevention of cross contamination by reducing lint

Resistant to liquids with its barrier feature

Resistant to moisture and bacterial penetration

Puncture, tear and abrasion resistant.

Free from toxic substances

Non-Sterile

Application:

This gown can be used during the surgery operations.

This product is disposable and it prevents the bacterial transfer between the medical team and the patient, reduces the occurring of possible contamination and infection between the patient and the medical team.

Precautions:

The product should be stored in a clean, dry place with proper ventilation and out of direct sunlight at 0°C-30°C (its relative temperature should not exceed 75%).